

How to pair your iPhone compatible hearing devices:

Step 1

- Go to **Settings** and select **General**

Step 2

- Select **Accessibility**

Step 3

- Select **Hearing Aids**

Step 4

- Place your Hearing Aids close to the iPhone, iPad or iPod.
- Open and close the battery door.
- The Apple device will then search for your hearing aids. Once connected, select the hearing aids under **Device**.

Step 5

- Select **Pair** for each hearing device. Wait 2 minutes before leaving screen.

Step 6

- Once complete the hearing aids will appear.
Congratulations! Your aids are now paired.